

Leksell GammaPlan®

- Share data easily
- Plan with flexibility
- Expand connectivity

# DICOM-RT

Experience increased versatility


# Experience increased versatility and smarter workflow

DICOM-RT is one of the many value-add extensions – add-ons – for Leksell GammaPlan®.


DICOM-RT add-on provides you with the opportunity to benefit from data generated by various treatment planning systems. It improves workflow as well as flexibility and facilitates highest possible clinical outcome by using all relevant treatment data – regardless type of treatment planning system.

Based on the DICOM-RT standard, the add-on provides you with significant advantages:

- Import user-defined volumes created on another treatment planning system. The volumes can be targets, organ at risk, isodoses or any other region of interest defined by the user.



Leksell GammaPlan doseplan




Leksell GammaPlan doseplan exported to Ergo evaluate planning software

- Export data created in Leksell GammaPlan. Share planning information such as outlined volumes and dose matrices with any DICOM-RT compatible system.
- Contour targets and organs at risk days in advance. The images and outlines are then imported through the DICOM-RT add-on into Leksell GammaPlan for dose planning.

Please contact your local Elekta representative for a hands-on demonstration of this versatile tool.

For further information go to [www.elekta.com/dicom-rt](http://www.elekta.com/dicom-rt).



## More power to your planning

Leksell GammaPlan functionality is continuously expanding. The latest software releases for Gamma Knife surgery are using state-of-the-art tools based on Linux and PC. This opens up for future upgrades and opportunities with each add-on meeting specific needs while giving improved workflow, cost-efficiency and treatment quality.

[www.elekta.com](http://www.elekta.com)

Human Care Makes the Future Possible

Corporate Head Office:

Elekta AB (publ)  
Box 7593, SE-103 93 Stockholm, Sweden  
Tel +46 8 587 254 00  
Fax +46 8 587 255 00  
[info@elekta.com](mailto:info@elekta.com)

Regional Sales, Marketing and Service:

North America,  
Atlanta, USA  
Tel +1 770 300 9725  
Fax +1 770 448 6338  
[info.america@elekta.com](mailto:info.america@elekta.com)

Europe, Latin America,  
Africa, Middle East & India  
Tel +44 1293 654068  
Fax +44 1293 654655  
[info.europe@elekta.com](mailto:info.europe@elekta.com)

Asia Pacific  
Hong Kong, SAR China  
Tel +852 2891 2208  
Fax +852 2575 7133  
[info.asia@elekta.com](mailto:info.asia@elekta.com)

