

One-Time Data Imports

DESCRIPTION

IMPAC's one-time data imports save valuable system configuration and start-up time by automatically uploading previous data captured from billing, practice management, EMR, hospital information systems or ad-hoc databases into the MOSAIQ system.

STANDARD IMPORTED DATA

The following standard one-time data imports are available at established rates and can be performed for new customer installations. Specifications for each import are available upon request.

IMPORTS	DESCRIPTION	FORMAT	COMMENTS
Payer	Insurance company ID, plan ID, address, contact information	Comma delimited	
Physician	Staff physicians and referring physician name, address, contact information, IDs, and licenses	Comma delimited	
Patient	Patient address, demographic data, next of kin, guarantor, insurance plans, emergency contact information, and diagnoses	Comma delimited	
Documents (Word, RTF)	Progress notes, nurses notes, prescriptions, Radiology reports, Pathology reports, and letters	Comma delimited, MS ACCESS	RTF Documents must have format characters that are supported by the IMPAC document viewer. To check the formatting, submit samples to your data import representative. Pre-requisites: Primary Patient Identifier shall exist in the IMPAC application (Patient Import, ADT interface or manual entry)
Scanned documents (tif, gif, bmp)	Insurance cards, referral slips, lab results, and images	Comma delimited, MS ACCESS	Pre-requisites: <ul style="list-style-type: none"> Primary Patient Identifier shall exist in the IMPAC application (Patient Import, ADT interface or manual entry)

Schedules	Patient appointments: Date, activity, location and staff	Comma delimited, MS ACCESS	Pre-requisites: <ul style="list-style-type: none"> Primary Patient Identifier shall exist in the IMPAC application (Patient Import, ADT interface or manual entry)
Patient Notes	Patient clinical notes: Clinical-family history, Clinical-oncology history, Clinical-past prescriptions, Clinical-social history, Clinical- Radiotherapy Summary	Comma delimited, MS ACCESS	Clinical notes will appear under eCHART Navigator Pre-requisites: <ul style="list-style-type: none"> Primary Patient Identifier shall exist in the IMPAC application (Patient Import, ADT interface or manual entry)
Diagnosis Detail	Diagnosis and Staging : Histology, Stage, TNM, Nodes	Comma delimited, MS ACCESS	Pre-requisites: <ul style="list-style-type: none"> Primary Patient Identifier shall exist in the IMPAC application. (Patient Import, ADT interface or manual entry)
Clinical Trials	Trial category, status, phase, disease site, restrictions, sponsor, patient enrolled date	Comma delimited, MS ACCESS	Pre-requisites: <ul style="list-style-type: none"> Primary Patient Identifier shall exist in the IMPAC application. (Patient Import, ADT interface or manual entry)
Treatment Summary	Prescription and dose site summary	Comma delimited, MS ACCESS	Pre-requisites: <ul style="list-style-type: none"> Primary Patient Identifier shall exist in the IMPAC application. (Patient Import, ADT interface or manual entry)
Treatment Detail	Planning information such as prescription, field definitions, setup information, images as well as treatment delivery data	SQL SERVER DB	Pre-requisites: <ul style="list-style-type: none"> Primary Patient Identifier shall exist in the IMPAC application. (Patient Import, ADT interface or manual entry) Images to be provided in DICOM format DICOM license (image import)

THE IMPORT PROCESS

1. The data to be imported shall be extracted by the customer or their representative and sent to IMPAC in ASCII character-delimited format or as a MS ACCESS or SQL Server DB according to IMPAC specifications.
2. IMPAC will import the data into your new MOSAIQ database. The data will first be imported to a practice database to allow inspection by your department. Once IMPAC receives verification and approval that the data has been properly imported into the practice database, it will be imported to the production clinical database.

CUSTOMER RESPONSIBILITY

It is the customer's responsibility to extract the data from the non-IMPAC system(s) or secure third party services to write scripts or reports to place the data in the IMPAC-specified format. In special cases IMPAC can extract data from standard non-proprietary database structures such as Microsoft Excel and Access. IMPAC will perform data analysis and discuss the analysis with the customer or their agent in order to resolve data ambiguity, perform any needed data scrubbing, and address issues with the supplied data.

Customer shall be responsible for participating in the validation process and correcting any data inconsistencies that may arise during the evaluation.

CUSTOM DATA IMPORT

For data not included in the standard one-time data imports, IMPAC can perform analysis and estimate the price to import data into specified IMPAC product data fields as a custom job.

DISCLAIMER

IMPAC can only guarantee full import of data delineated in its import specifications. The specifications are available for review upon request.

IMPAC Medical Systems, Inc., IMPAC, the IMPAC logo, and our product names are trademarks of IMPAC. All other brand names or trademarks are the property of their respective holders. © 2009 IMPAC Medical Systems, Inc.