


Create Clinical Decision Support Reference

Create Clinical Decision Support (CDS) reference to be presented when a CDS Alert is generated

- 1) Click **File > Libraries > CDS Reference**
- 2) Click **Add**



The screenshot shows the 'CDS References Management' window. A table lists various CDS references with columns for Title, Reference, and Created By. A red box with the number '2' highlights the 'Add' button in the top-left corner of the window. Overlaid on the table is a 'CDS Reference' dialog box. A red box with the number '3' highlights the 'Caption' field in the dialog box, which contains the text 'COVID-19 Screening Protocol'. Below the caption field, there are sections for 'Therapeutic rationale:', 'Links:', and 'Reference:'. The 'Links:' section includes a list of items: 'a. Bibliographic citation:', 'b. Developer of the intervention:', 'c. Funding source of the intervention:', and 'd. Release:'. Below these is the 'Effective date:' field. At the bottom of the dialog box, there is a red text label that reads 'CDS Reference will be added'.

- 3) Enter **Caption**
- 4) Optionally, enter instructions and bibliographic reference with or without URL links.
- 5) Click **OK** to save