



MOSAIQ IQ Script: Import a Screening Assessment

Import a COVID-19 screening Assessment form


- 1) Download Assessment **Optional: Download Custom Layout**
- 2) Save to an accessible location
- 3) In MOSAIQ, select **File > Import/Export > Observation View > Import**
- 4) To import Assessment: Select **Open File**, navigate to **COVID19_Assessment_<date>.xml**



- 5) Reconcile any conflicts




- 6) Select checkbox, select **Import**



- 7) Add or remove content using **Chart Builder** (see Quick Guides)
- 8) Configure a **Merge Field** or **Merge Field Widget** if desired (see Quick Guides)

Optional: Apply a Custom Layout


- 1) Double-click COVID-19 Assessment, select Layout



MOSAIQ IQ Script: Import a Screening Assessment continued

Import a COVID-19 screening Assessment form

- 2) Select **Layout Options > Customize Layout**. Search folder and navigate to **COVID19_Layout_<date>.xml**. Resize if desired and **Save**.



Custom Layout:

The screenshot shows the 'COVID-19 Screening Assessment' form in the 'Custom Layout' view. The form is titled 'COVID-19 Screening Assessment' and includes a 'Record Date and Time' field showing '5/4/2020 6:36:34 AM'. The form is divided into several sections:

- Symptoms:** A grid of checkboxes for various symptoms. A red box with the number '14' is overlaid on the 'None Reported' checkbox. The symptoms listed are: None Reported, Cough, Sore Throat, Loss of Smell or Taste, Fever, Chills or Sweating, Difficulty Breathing, Body Aches, Vomiting or Diarrhea, and Headache.
- Risk Factors:** A series of dropdown menus for 'Exposed to COVID-19?', 'Medical, Other High Risk Work?', 'Travel in Last 14 Days', and 'Travel Location(s)'. There is also a text input field for 'Comment - High Risk Work'.
- Precautionary Measures:** A series of dropdown menus for 'Practicing Social Distancing?', 'Mask in Public', 'Isolating Household', and 'Reduced use of public spaces?'. There is also a text input field for 'Other Precautionary Measures'.
- Add Additional Note:** A text input field for 'Comment - COVID-19:'.

At the bottom of the form, there is a red status bar that reads 'COVID-19 Screening Assessment Will Be Added'.